

# G-FORCE Braces

## Cryo-Pneumatic Knee Brace

Latex Free

### Directions for use of the Cryo-Pneumatic Knee Brace:

1. Remove optional gel pack and place in freezer for a minimum of two hours. **(PICTURE 1A)**
2. Remove gel from freezer, re-attach gel pack inside knee brace, making sure seams line up over kneecap (patella) opening. **(PICTURE 1B)**
3. Grasp knee brace and locate top (valve located at top) and place opening over affected kneecap (patella). **(PICTURE 1C)**
4. Secure 2 compression straps to desired fit. **(PICTURE 2)**
5. Attach tube with bulb valve. **(PICTURE 3A)**
6. To inflate, switch to "ON" position, pump bulb to desired compression, taking care not to overinflate. To deflate, press button at base of bulb.
7. When desired compression is attained, switch to "OFF" position and remove tube with bulb from valve connection.

Note: brace may be deflated by removing tube while in "ON" position.


**WARNING: Recommended cold therapy treatment: 20 minutes on patient, 20 minutes off**



### CAUTION:

1. **DO NOT OVER INFLATE** – improper inflation may inhibit circulation or cause skin irritation
2. Use only after following directions for use
3. Do not use gel pack on area of sensitive skin
4. Cold therapy should not be used by people with circulatory problems unless prescribed by a physician
5. Gel pack is non-toxic. Not intended for ingestion
6. Discard gel pack if punctured or torn

**Contents: Device**



PICTURE 1



PICTURE 2



PICTURE 3

**G-FORCE Braces**

300 Bedford Street, Manchester, New Hampshire 03101

[www.gforcebraces.com](http://www.gforcebraces.com)

Toll Free (888) 248-9924