

PRODUCT CATALOG



www.elitemedicalsupply.com

ELITE  MEDICAL[®]
SUPPLY



We provide high-quality orthopedic braces and bone growth stimulators to patients being treated for back, knee, and extremity conditions.

P 866-712-0881

F 866-712-0882

To Order: order@elitemedicalsupply.com

Questions: info@elitemedicalsupply.com

1900 Ridge Road, Suite 125

West Seneca, NY 14224



WHO WE ARE

Elite Medical Supply was founded over 20 years ago just outside Buffalo, NY. We have specialized in back, knee, neck, and extremity pain disorders since 1998. Our staff of Chiropractors and Certified Orthotic Fitters are experts at identifying the equipment to best serve your specific needs.

We are in-network with over 150 insurance plans, Medicare, Medicaid, Workers' Compensation and No-Fault. Click [here](#) to view our in-network insurance plans.

Whether for post-operative rehabilitation, recovery following an accident, or in cases of chronic pain, our staff will work collaboratively with your healthcare provider to find the correct pain management solution.

We are a Medicare-Enrolled Supplier of DMEPOS

Elite Medical Supply is a Medicare-enrolled supplier of durable medical equipment that accepts Medicare assignment and bills Medicare. We have grown from servicing New York State beneficiaries to participating in nationwide insurance networks.

Click below to learn more about Elite Medical Supply





HOW TO ORDER

How simple can we make the ordering process?

What if all you had to do was send a prescription and supporting documents?

Well that's how easy it is to order with us!

Use our online form, fax, or email a *prescription and supporting documents.

www.elitemedicalsupply.com/referrals

866-712-0882

order@elitemedicalsupply.com

To expedite order processing, our team needs relevant patient information:

- Medical records supporting necessity of product ordered
- Demographic Sheet
- Insurance Information
- Waist measurement around naval for spinal products
- Measurements 6" above and below knee joint space for knee products

*Include the diagnosis code(s) on the prescription

Elite Medical Supply

1900 Ridge Road

Suite 125

West Seneca, NY 14224

866-712-0881



PRODUCT LINKS

BACK BRACES

KNEE BRACES

LOWER EXTREMITY

UPPER EXTREMITY

BONE GROWTH STIMULATORS

Visit www.elitemedicalsupply.com to view all of our products.



BACK BRACES



Our back brace line includes bracing products that cover the entire spine.

LSO - Lumbar Sacral Orthosis

TLSO - Thoracic Lumbar Sacral Orthosis

SI - Sacroiliac Belt

Posture Correction

BACK BRACES

Elite Air Decompression LSO

The Elite Air Decompression LSO (lumbar sacral orthosis) combines the support of rigid panels with a decompression lumbar belt to support and relieve pain. The front and back panels are removeable to offer continued therapy options during recovery. The back panel extends from the sacrococcygeal junction to the T-9 vertebra. Decompression of the lower spine is accomplished through expansion of the lumbar belt from 5.1" to 7.9" (13 – 20cm) which aids in relaxing muscle spasms and relieves pain.

AIR DECOMPRESSION LSO WITH RIGID SUPPORT PANELS & STEP-DOWN SUPPORT



KEY FEATURES

- Adjustable air decompression
- Rigid front and back (anterior and posterior) panels for superior support
- Removeable panels for step-down support
- Prevents and relieves back pain, relaxes muscle spasms
- Belt expands from 5.1" to 7.9"
- Sizes from 26" to 58" (Small - 3XLarge)
- Extender panel increases the circumference up to 8"
- HCPCS L0648

TYPICAL INDICATIONS

- Burst fracture
- Bulging or herniated discs
- Degenerative disc disease
- Post-operative recovery
- Sciatica
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis

BACK BRACES

Aspen OTS™ LSO 650

The Aspen OTS™ LSO 650 is a one-size adjustable, inelastic spinal brace featuring removable lateral panels that provide increased control for patients in need of additional support. This clinically-tested, non-narcotic, non-invasive treatment solution combines support and direct independent compression to promote effective trunk stability, enabling patients to resume activities of daily living. Comfortably fits waist sizes ranging from 28-54 inches, and up to 64 inches with Universal Extension Panel.

ONE SIZE ADJUSTABLE INELASTIC SPINAL BRACE



KEY FEATURES

- Low profile design
- One size adjustable
- Targeted compression system
- Firm and flexible back panel
- HCPCS L0650

TYPICAL INDICATIONS

- Chronic Back Pain
- Degenerative, Herniated or Bulging Disc
- Spondylolysis, Spondylolisthesis
- Spinal Stenosis
- Facet Syndrome
- Tumors
- Spinal Instability, Burst FX
- Pre or Post-Op Spinal Lumbar Stabilization
- Pre or Post-Op Disc Surgery

BACK BRACES

DJO EXOS® FORM™ II 637

The Exos FORM™ II 637 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos FORM™ II 637 provides superior sagittal and coronal control and support from T9-S1. The BOA® Fit System is a dial-based performance fit system that is engineered with high quality, durable materials that enable a microadjustable connection that's built to perform.

DIAL IN TO FAST, EFFORTLESS, PRECISION FIT



KEY FEATURES

- V-STX™ Material
- ComfortCORE™ Foam
- Diamond2 Grid™
- BOA® Fit System - fast, effortless, precision fit
- Adjustable Belt Wings
- Modular Design with step-up/step-down design
- Semi-Universal Sizing
- Also available as HCPCS L0631, L0627, L0626, L0621

TYPICAL INDICATIONS

- Sacroiliac Pain
- Acute and Chronic Lower Back Pain
- Lower Back Sprains / Strains
- Lumbar Disc Displacement
- Osteoporosis
- Disc Herniation and Degeneration
- Post-Operative Laminectomy
- Post-Operative Disectomy
- Spondylolisthesis/Spinal Stenosis

BACK BRACES

Comfortland Delta 37

Comfortland's Delta 37 is the perfect combination of comfort and support, giving practitioners a cost-effective LSO without compromising construction and functionality. Weighing in at 1.5lbs, the Delta 37 puts little strain on the wearer, while a revolutionary pulley compression system delivers a 4:1 return on the user's pulling power. Universal sizing accommodates waist sizes from 25" to 68", and a sturdy 15" posterior panel and lateral support panels help prevent any abnormal movements of the lower spine.

PERFECT COMBINATION OF COMFORT & SUPPORT



KEY FEATURES

- Universal size
- Optional 8" extender panel available
- 4:1 Compression ratio
- Heat moldable posterior panel
- Removable lateral panels
- HCPCS L0637/L0650

TYPICAL INDICATIONS

- Spinal instability
- Spinal stenosis
- Chronic back pain

BACK BRACES

Aspen OTS™ LSO 648

The Aspen OTS™ LSO 648 is an adjustable, inelastic spinal brace offering a clinically-tested technology option for patients experiencing low back muscle spasms causing pain. This non-narcotic, non-invasive solution combines support and direct independent compression to promote effective trunk stability, enabling patients to regain activities of daily living. Comfortably fits waist sizes ranging from 28-54 inches, and up to 64 inches with Universal Extension Panel.

ONE SIZE ADJUSTABLE INELASTIC SPINAL BRACE



KEY FEATURES

- Low profile design
- One size adjustable
- Targeted compression system
- Firm and flexible back panel
- HCPCS L0648

TYPICAL INDICATIONS

- Chronic Back Pain
- Degenerative, Herniated or Bulging Disc
- Spondylolysis, Spondylolisthesis
- Spinal Stenosis
- Facet Syndrome
- Tumors
- Spinal Instability, Burst FX
- Pre or Post-Op Spinal Lumbar Stabilization
- Pre or Post-Op Disc Surgery

BACK BRACES

Optec Venum Hybrid LSO & TLSO

The Venum Hybrid LSO is a one-size brace that simplifies fitting while providing the support and stabilization of a traditional lumbar brace. This one-of-a-kind brace also provides unique support through a dual drawstring pulley closure to increase the support and compression delivered to the low back with minimal effort to the spine from the T-9 vertebra to the sacrococcygeal junction. The detachable polymer panels allow for easy adjustment of the compression level and support throughout rehabilitation.

ONE-SIZE LUMBAR BRACES WITH ADAPTABLE STEP-DOWN SUPPORT



KEY FEATURES

- Rigid Panels and Step-Down Support
- Made from lightweight and breathable mesh
- Detachable polymer panels allow for easy adjustment of compression and support
- One-size brace simplifies fitting
- Dual drawstrings increase support and compression

TYPICAL INDICATIONS

- Burst Fracture
- Bulging or Herniated Discs
- Degenerative Disc Disease
- Post-Op
- Sciatica
- Spinal Stenosis
- Spondylolysis
- Spondylolisthesis

BACK BRACES

Venum Hybrid Pro

The Venum Hybrid Lite LSO is the ultimate one-size-fits-all soft comfort brace for the lumbar-sacral region. Its primary function is to reduce pain, provide support and immobilization to the lumbar-sacral spine. The Venum Hybrid design allows each patient to get the perfect compression and immobilization for their recovery. The built-in anterior panels produce intra-abdominal pressure to reduce the load on discs thus reducing pain and improving post-operative outcome when used appropriately.

ONE SIZE FITS ALL COMFORT BRACE



KEY FEATURES

- Adjustable One-Size-Fits-All Design
- Step-down, streamlined design
- Dual compression closure system provides 6:1 pull ratio
- Ergonomic arthritic grips
- Made of durable, highly breathable, soft, hypoallergenic materials that wick away perspiration
- Fits waist circumferences from 21" – 63"
- Optional extender panel adds 10" to circumference
- HCPCS L0631/L0648

TYPICAL INDICATIONS

- Chronic Back Pain
- Degenerative or herniated or bulging disc
- spondylolysis, spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumors
- Spinal instability, burst FX
- Pre or post-op spinal lumbar stabilization
- Pre or post-op disc surgery

BACK BRACES

Ovation Medical Nu-Form Universal Back Brace

The Nu-Form universal back brace has unmatched breathability and durability. The Nu-Form back brace design includes an innovative pressure relief zone and a uniquely tailored fit. Featuring Fusion Form technology that results in a stronger, lighter, and more breathable brace.

INNOVATIVE PRESSURE RELIEF ZONE WITH LOW-PROFILE DESIGN



KEY FEATURES

- Innovative pressure relief zone
- Ultra-low profile
- Universal sizing
- 5 year warranty
- HCPCS L0637 or L0650
- Also available as HCPCS L0642, L0648

TYPICAL INDICATIONS

- Burst Fx
- Bulging or Herniated Discs
- Degenerative Disc Disease
- Post-Op
- Sciatica
- Spinal Stenosis
- Spondylolysis
- Spondylolisthesis

BACK BRACES

Comfortland Delta 31

Comfortland's Delta 31 is the perfect combination of comfort and support, giving practitioners a cost-effective LSO without compromising construction and functionality. Weighing in at 1.5lbs, the Delta 31 puts little strain on the wearer, while a revolutionary pulley compression system delivers a 4:1 return on the user's pulling power. Universal sizing accommodates waist sizes from 25" to 68", and a sturdy 15" posterior panel helps prevent any abnormal movements of the lower spine.

REVOLUTIONARY PULLEY COMPRESSION SYSTEM



KEY FEATURES

- Universal Size
- Optional 8" extender panel available
- Breathable mesh construction
- 4:1 Compression ratio
- Heat moldable posterior panel
- HCPCS L0631/ L0648

TYPICAL INDICATIONS

- Spinal instability
- Spinal stenosis
- Chronic back pain

BACK BRACES

Aspen OTS™ Lumbar 642

The Aspen OTS™ Lumbar 642 features the lowest profile design within the Aspen OTS line. This one-size adjustable, inelastic brace offers a clinically-tested technology for patients experiencing low back muscle spasms causing pain. As a non-narcotic, non-invasive treatment solution, it combines the needed level of support coupled with direct independent compression to promote effective trunk stability, enabling patients to regain activities of daily living. Comfortably fits waist sizes ranging from 28-54 inches.

LOW-PROFILE, ONE-SIZE ADJUSTABLE



KEY FEATURES

- Low profile design
- One size adjustable
- Fits up to 64 inches with Universal Extension Panel
- Targeted compression system
- HCPCS L0642

TYPICAL INDICATIONS

- Chronic Back Pain
- Degenerative, Herniated or Bulging Disc
- Spondylolysis, Spondylolisthesis
- Spinal Stenosis
- Facet Syndrome
- Tumors
- Spinal Instability, Burst FX
- Pre or Post-Op Spinal Lumbar Stabilization
- Pre or Post-Op Disc Surgery

BACK BRACES

Optec Venum Hybrid Lite LSO

The Venum Hybrid Lite LSO is the ultimate one-size-fits-all soft comfort brace for the lumbar-sacral region. Its primary function is to reduce pain, provide support and immobilization to the lumbar-sacral spine. The Venum Hybrid design allows each patient to get the perfect compression and immobilization for their recovery. The built-in anterior panels produce intra-abdominal pressure to reduce the load on discs thus reducing pain and improving post-operative outcome when used appropriately.

ONE SIZE FITS ALL LSO WITH PERFECT COMPRESSION



KEY FEATURES

- Adjustable One-Size-Fits-All Design
- Step-down, streamlined design
- Dual compression closure system provides 6:1 pull ratio
- Ergonomic arthritic grips
- Made of durable, highly breathable, soft, hypoallergenic materials that wick away perspiration
- Fits waist circumferences from 21" – 63"
- Optional extender panel adds 10" to circumference
- HCPCS L0627/L0642

TYPICAL INDICATIONS

- Chronic Back Pain
- Degenerative or herniated or bulging disc
- Spondylolysis, spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumors
- Spinal instability, burst FX
- Pre or post-op spinal lumbar stabilization
- Pre or post-op disc surgery

BACK BRACES

Comfortland Delta 27

Comfortland's Delta 27 is the perfect combination of comfort and support, giving practitioners a cost-effective LSO without compromising construction and functionality. Weighing in at 1.5lbs, the Delta 27 puts little strain on the wearer, while a revolutionary pulley compression system delivers a 4:1 return on the user's pulling power. Universal sizing accommodates waist sizes from 25" to 68", and a sturdy 10" posterior panel helps prevent any abnormal movements of the lower spine.

REVOLUTIONARY PULLEY SYSTEM



KEY FEATURES

- Universal size
- Optional 8" extender panel available
- Breathable mesh construction
- 4:1 Compression ratio
- Heat moldable posterior panel
- HCPCS L0627/ L0642

TYPICAL INDICATIONS

- Spinal instability
- Spinal stenosis
- Chronic back pain

BACK BRACES

Comfortland Avalon LSO

The Comfortland Avalon LSO is a great economic solution to low back support. Comprised from a breathable mesh fabric and featuring a 14" posterior panel with "classic" elastic compression pulls, the Avalon LSO is a great lower spine support designed to fit in any budget.

ECONOMIC LSO WITH CLASSIC FEATURES



KEY FEATURES

- Low profile belt design
- Lightweight and breathable construction
- Removable 14" posterior panel
- HCPCS L0643

TYPICAL INDICATIONS

- Spinal instability
- Spinal stenosis
- Chronic back pain

BACK BRACES

Comfortland Delta 26

Comfortland's Delta 26 is the perfect combination of comfort and support, giving practitioners a cost-effective back brace without compromising construction and functionality. The Delta 26 puts little strain on the wearer, while universal sizing accommodates waist sizes from 25" to 68", and a sturdy 15" posterior panel helps prevent any abnormal movements of the lower spine.

PERFECT COMBINATION OF COMFORT & SUPPORT



KEY FEATURES

- Universal size
- Optional 8" extender panel available
- Breathable mesh construction
- Pulley system to maximize compression
- Easy to use pull tabs
- HCPCS L0626/L0641

TYPICAL INDICATIONS

- Spinal instability
- Spinal stenosis
- Chronic back pain

BACK BRACES

Aspen OTS™ TLSO 457

The Aspen OTS™ TLSO 457 is a one-size adjustable, inelastic spinal brace designed to provide additional thoracolumbar support, bringing patients back to a functional midrange to help create a dynamic environment for healing from complex spinal conditions. This clinically-tested, non-narcotic, non-invasive treatment solution combines support and direct independent compression to promote effective trunk stability, enabling patients to regain activities of daily living.

ONE-SIZE ADJUSTABLE, INELASTIC SPINAL BRACE



KEY FEATURES

- Designed for extended wear
- Low profile
- Pull-through sizing mechanism
- Aspen's patented 4:1 tightening mechanism
- An optional Extension Panel is available that increases circumference by 10"
- HCPCS L0457

TYPICAL INDICATIONS

- Postoperative support / pain relief
- Kyphosis secondary to osteoporosis
- Decompressive procedures
- Degenerative disc disease
- Bulging or herniated disc
- Fracture management
- Structural kyphosis
- Spondylolisthesis/Spondylosis
- Spinal stenosis

BACK BRACES

Aspen Active™ P-TLSO

The Aspen Active™ Postural-TLSO is a lightweight, adjustable brace designed to address the symptoms associated with poor posture often due to prolonged use of electronic devices/technology. While actively influencing a patient's posture, the brace provides gentle support while activating postural stabilizing muscles. This product helps train the body to hold neutral posture by concurrently stretching and strengthening the structures that have become tight and deconditioned over time due to poor posture.

THE NEWMATIX™ SYSTEM ACTIVATES MUSCLES TO ACHIEVE A NEUTRAL POSITION



KEY FEATURES

- NEWmatix™ System activates muscles to achieve a neutral position
- One size adjustable
- Low-profile
- Shoulder straps are adjustable to fit various bodies
- Independent tightening mechanism

TYPICAL INDICATIONS

- Postoperative support / pain relief
- Kyphosis secondary to osteoporosis
- Decompressive procedures
- Degenerative disc disease
- Bulging or herniated disc
- Fracture management
- Structural kyphosis
- Spondylolisthesis/Spondylosis
- Spinal stenosis

BACK BRACES

Comfortland Delta TLSO with ATE

The Delta TLSO with ATE (anterior thoracic extension) is designed to provide motion restriction for the thoracic and lumbar regions. The posterior extends from the sacrococcygeal junction to just inferior of the scapular spine with tongues providing lateral strength. The Anterior extends from the symphysis pubis to the sternal notch to restrict forward flexion and encourage proper posture. This user-friendly TLSO allows you to restrict the thoracic vertebra through precise compression with its unique closure system.

THORACIC & LUMBAR REGION MOTION RESTRICTION



KEY FEATURES

- Universal size
- Adjustable anterior bar
- Adjustable thoracic support
- Lightweight, durable construction
- Optional 8" extender panel available
- HCPCS L0456/L0457

TYPICAL INDICATIONS

- Sprains and Strains of the thoracic region
- Displacement of intervertebral discs,
- Fracture of the vertebral column
- Burst fractures
- Herniated or bulging discs
- Post-op spinal stabilization
- Spinal stenosis

BACK BRACES

Comfortland Delta 12X

Designed to aid with compression fractures, upper back pain, osteoporosis, and kyphosis, the Delta 12X is the first TLSO in the Delta series. With an adjustable posterior support, adjustable shoulder straps, and the unique abdominal compression system native to the Delta braces, the DL-12X helps reduce daily pain, increase core strength, and improve the quality of life for those who suffer from chronic thoracic pain.

ADJUSTABLE POSTERIOR SUPPORT AND SHOULDER STRAPS



KEY FEATURES

- Universal size
- Optional 8" extender panel is available
- 4:1 Compression ratio
- Adjustable thoracic support
- Lightweight and durable construction
- HCPCS L0456/L0457

TYPICAL INDICATIONS

- Sprains and Strains of the thoracic region
- Displacement of intervertebral discs
- Fracture of the vertebral column

BACK BRACES

Comfortland Delta 18X TLSO

Designed for back pain, sprains and strains of the cervical and thoracic regions, and for prophylactic use to prevent injury, the Delta 18X TLSO is crafted from high-quality soft goods, giving optimum compression and support while remaining lightweight and comfortable. Providing a highly adjustable, universal support, the Delta 18X fits most patients with ease, drastically reducing ordering costs and stock.

OPTIMUM COMPRESSION & SUPPORT



KEY FEATURES

- High quality soft goods giving optimum compression
- Light-weight and comfortable
- Adjustable and universal support
- HCPCS L0450

TYPICAL INDICATIONS

- Sprains and strains of the cervical and / or thoracic region
- Prophylactic use to prevent injury
- Support for weak muscles
- Spinal stenosis of thoracic region

BACK BRACES

Comfortland Delta 455X TLSO

Indicated for sprains and strains of the thoracic region, back pain, and for prophylactic use, the Comfortland Delta 455X TLSO is designed to provide motion restriction for the thoracic and lumbar regions. The posterior extends from the sacrococcygeal junction to just inferior of the scapular spine . The anterior extends from the symphysis pubis to the sternal notch to restrict forward flexion and encourage proper posture. This user-friendly TLSO allows you to restrict the thoracic vertebra through precise compression with its unique closure system.

MOTION RESTRICTION FOR THE THORACIC AND LUMBAR REGIONS



KEY FEATURES

- Universal sizing
- Adjustable anterior bar
- Adjustable thoracic support
- Light-weight and durable construction
- HCPCS L0455

TYPICAL INDICATIONS

- Sprains and strains of the thoracic region
- Back pain
- Prophylactic use

BACK BRACES

Aspen Evergreen™ 621 SI Belt

The Evergreen™ 621 SI belt utilizes a bi-lateral pull system that provides maximum compression without the rotational loading that can occur with a single pull tab system. The low profile and lightweight design makes extended wear comfortable while maintaining the quality of design and durability that you've come to expect from Aspen.

MAXIMUM COMPRESSION WITHOUT ROTATIONAL LOADING



KEY FEATURES

- Low Profile and Lightweight
- Easy to pull regardless of patient strength
- Silicone strips on the inside of the brace hold fast and prevent migration especially while moving and sitting
- Utilizes Aspen's proven SlickTrack™ mechanically advantaged tightening system
- Durable two-ply construction with breathable mesh fabric for comfortable and extended wear

TYPICAL INDICATIONS

- Fracture management
- Spinal stenosis
- Sprains / strains
- Facet syndrome
- SI dysfunction
- Radiculopathy
- Spondylosis
- Sacroiliitis
- Injections

Comfortland Sacro Stabilizer

The Comfortland SI-23 Sacro Stabilizer is designed to provide support and relief to abdomen and sacroiliac joints. Featuring an easy to adjust universal sizing system and both anterior and posterior panels, the SI-23 is a great solution for SI issues.

ABDOMEN & SACROILIAC JOINT SUPPORT



KEY FEATURES

- Customizable universal design
- Easy to use hook and loop fastening
- Soft, breathable construction
- Rigid anterior and posterior panels provide support and relief to abdomen and sacroiliac joints
- HCPCS L0623

TYPICAL INDICATIONS

- SI joint instabilities
- SI joint arthritis
- SI joint syndrome
- SI joint pain

BACK BRACES

Comfortland SI 100 Belt

The Comfortland SI-100 is designed to provide immediate pain relief and quickly restore functional mobility. With no irritating buckles or pads, the SI-100 utilizes elastic bands to provide compression and can be worn for extended periods of time without causing muscle weakness or atrophy.

IMMEDIATE PAIN RELIEF WITHOUT MUSCLE WEAKNESS OR ATROPHY



KEY FEATURES

- Low-profile design
- Adjustable Velcro enclosures
- Designed to be worn for extended periods without causing muscle atrophy
- HCPCS L0621

TYPICAL INDICATIONS

- Piriformis syndrome
- Hip/pelvic/upper leg pain
- Low back pain

KNEE BRACES



Our knee brace line includes bracing products for various conditions.

ACL - Anterior Cruciate Ligament injuries

OA - Osteoarthritis conditions

ROM - Range of Motion

Hinged

Immobilizer

Knee Cage

Post-Op

Patella Stabilizer

KNEE BRACES

Aspen OA Knee+

This premium knee brace is ideal for patients with mild to moderate osteoarthritis. The low profile, universal design allows for less inventory and better patient compliance. Adjustable stops are included to set the range of motion for appropriate treatment during continuum of care.

UNIVERSAL FIT AND LOW PROFILE



KEY FEATURES

- Two-point hinge adjustment
- Range of motion adjustable stops
- One-step closure
- Comfortable, non-slip pads help to prevent migration
- Low-profile design
- Left/right specific
- Universal Size
- HCPCS Code L1843, L1851

TYPICAL INDICATIONS

- For persons with unicompartement knee osteoarthritis (OA), off-unloader braces are mechanical intervention designed to reduce pain, improve physical function, and possibly slow disease progression.

KNEE BRACES

Ovation Medical Game Changer OA Knee Brace

The Game Changer OA knee brace offers superior relief for osteoarthritic pain. The single upright hinge design is comfortable, lightweight, and functional. Treats both medial and lateral knee pain while being easy to apply.

TREAT MEDIAL AND LATERAL KNEE PAIN



KEY FEATURES

- Superior offloading
- OA knee pain relief
- Low profile
- Highly functional
- Left/right specific
- Universal size

TYPICAL INDICATIONS

- For persons with unicompartement knee osteoarthritis (OA), off-unloader braces are mechanical intervention designed to reduce pain, improve physical function, and possibly slow disease progression.

KNEE BRACES

Comfortland OA Knee Brace

Designed for patients with mild to moderate osteoarthritis of the knee, the OA-200 is a low-profile brace with 3 points of leverage for reducing pressure inside the knee joint. It's universal; fitting a wide variety of patients easily. The strapping system allows for full adjustment, but also features quick snaps for quick and easy application.

UNIVERSAL FIT WITH QUICK SNAPS FOR EASY DONNING



KEY FEATURES

- Universal size
- Medial or lateral unloading for osteoarthritic conditions
- Breathable liners laminated with silicone eliminate migration
- Optimal brace for protection and support of ligament injuries and instabilities
- HCPCS L1843

TYPICAL INDICATIONS

- Provides stability and support following injury or reconstruction of the ACL, PCL, MCL, LCL and menisci.

KNEE BRACES

Aspen ROM Knee Brace

The Aspen ROM Knee is a range of motion brace that provides optimal compression and support. The brace comes in both sized and universal options. The sized version features a premium AirX material that provides a high degree of breathability and moisture wicking, allowing for extended wear and improved patient compliance. The Universal option provides adjustability to fit most patients. Made from a lightweight neoprene material, this option allows for optimal compression and comfort.

UNIVERSAL SIZE WITH OPTIMAL COMPRESSION & SUPPORT



KEY FEATURES

- Easy access range of motion hinge
- Posterior strap for additional adjustment and comfort
- Sturdy hinge struts for optimal customization
- Lightweight neoprene material for ideal compression and breathability
- Universal size
- HCPCS L1833

TYPICAL INDICATIONS

- Patellar Honodroopathy
- Retropatellar Osteoarthritis, post-operative and post traumatic irritations including rehabilitation of: MCL, LCL, PCL, ACL
- Post lateral release surgery to limit the extension and flexion of the joint
- Patella and knee instabilities
- General knee pain, Mild Chondronmalacia
- Condylar/Tibial Plateau Fractures

KNEE BRACES

Optec ROM Patella Stabilizer (ROMPS)

The ROMPS knee brace was built to offer durable, dependable support to patients of every activity level. With various heights and designs available, the ROMPS knee brace ensures that patients will receive the optimal brace for their specific needs. The ROMPS knee brace brings relief to patients with minor ligament issues, patella tracking issues and other patella related diagnoses with horse-shoe padding to keep the patella in place.

CONTROLS FLEXION AND EXTENSION WITH PATENTED, ADJUSTABLE ALUMINUM RODS



KEY FEATURES

- Horseshoe padding to keep the patella in place
- Stabilize & comfortably contain the patellar compartment
- Control flexion & extension with patented, adjustable aluminum rods
- Support for patients of every activity level

TYPICAL INDICATIONS

- Patellar Honodroopathy
- Retropatellar Osteoarthritis
- Post-op and post-traumatic irritations including rehab of: MCL, LCL, PCL, ACL
- Post lateral release surgery to limit the extension and flexion of the joint,
- Patella and knee instabilities
- General knee pain, Mild Chondronmalacia
- Condylar Fractures,Tibial Plateau Fractures

KNEE BRACES

Comfortland ACL Hinged Knee Brace

Designed to provide stability and support following injury or reconstruction of the ACL, PCL, MCL, and LCL, The CK-207 features a lightweight sturdy frame lined with durable soft goods for the combination of optimum support and comfort. The user friendly ROM hinge helps patients easily and quickly adjust the desired range of flexion and extension, making it the perfect brace for those active patients needing additional support during recovery.

LIGHTWEIGHT, STURDY FRAME FOR OPTIMUM SUPPORT



KEY FEATURES

- Lightweight, durable construction
- Malleable uprights
- Drop-lock ROM hinge
- HCPCS L1852

TYPICAL INDICATIONS

- Internal derangement of the knee
- Loose body in the knee
- Old disruption of anterior cruciate ligament
- Tear of Medial cartilage or meniscus of the knee
- Tear of lateral cartilage or meniscus of the knee
- Sprains and strains of the medial or lateral collateral ligament

KNEE BRACES

Comfortland Universal Hinged Knee Brace

The Comfortland Universal Hinged Knee Brace provides stability and support for chronic knee pains and instability. Available as a universal fit brace with a hidden hinge. Effective in treating mild injuries of the collateral ligaments, menisci, and patellar retinaculum.

UNIVERSAL KNEE BRACE WITH HIDDEN HINGE



KEY FEATURES

- Universal Size
- Lightweight Construction
- Breathable materials keep patients comfortable
- Helps restore functionality

TYPICAL INDICATIONS

- Effective in treating mild injuries of the collateral ligaments, menisci, and patellar retinaculum.
- Also internal derangement of the knee, chondromalacia of the patella, old disruption of the anterior cruciate ligament, chondromalacia, tear of the cartilage or meniscus of the knee, sprains and strains of the knee.

Comfortland Neoprene Hinged Knee Support

Indicated for mild or moderate ligament joint pain, pain and sprains in the knee, and Osgood Schlatter disease, the neoprene hinged knee is ideal for patients who prefer slip-on style knee supports. With a neoprene/4 way stretch fabric construction, easily removable hinges, as well as thigh and calf straps, the CK-105 is comfortable, easy to use, and provides optimum knee support.

EASY PULL-UP DESIGN WITH REMOVABLE HINGES



KEY FEATURES

- Easy pull-up design
- Sturdy lateral supports
- 4-way stretch nylon for patient comfort
- Easily removable hinges
- HCPCS L1812

TYPICAL INDICATIONS

- ACL, MCL, PCL, and LCL injuries
- Soft tissue injuries
- Post-operative knee immobilization

Optec Gladiator 1832 Knee Brace

The GLADIATOR 1832 Knee Brace can help speed recovery and is adjustable for a comfortable fit. It is recommended for noncontact sports and daily activities. The Gladiator 1832 provides compression and support and is ideal for speeding recovery of ACL and MCL injuries. The Knee brace is meant for patients who can load weight on the injured knee and is capable of some ambulation. The polycentric hinge has aluminum uprights and offers full flexion and extension control with settings at 0°, 15°, 30°, 60° and 90°.

UNIVERSAL OA KNEE BRACE WITH FULL STRAPPING ADJUSTMENT



KEY FEATURES

- Universal size
- O.M. hinges with contoured condyle shells
- Full flexion and extension control
- ROM settings at 0°, 15°, 30°, 60° and 90°
- HCPCS L1832, L1833

TYPICAL INDICATIONS

- Patellar Honodroopathy
- Retropatellar Osteoarthritis, post-operative and post traumatic irritations including rehabilitation of :MCL, LCL, PCL, ACL
- Post-lateral release surgery to limit the extension and flexion of the joint, Patella and knee instabilities
- General knee pain, Mild Chondronmalacia
- Condylar Fractures, Tibial Plateau Fractures

KNEE BRACES

Breg T-Scope® Premier Knee Immobilizer

The Breg T Scope® Premier Post-Op knee brace provides support and protects range-of-motion for patients who have undergone knee surgery, sustained a knee injury, or are experiencing instability. A patented hinge controls the range of motion of the patient's knee with an easy-to-use extension lock that allows the brace to operate in five typical positions of knee extension. This unique knee brace has telescoping calf and thigh sleeves for adjustable sizing, plus trimmable straps and lock clips so patients can customize and secure the right fit.

ROM EXTENSION LOCK AND UNIVERSAL FIT



KEY FEATURES

- Easy to use ROM hinge with quick adjusting flexion and extension stops
- Telescoping calf and thigh sleeves for sizing a wide range of patients from 5' to 6' 4" tall
- Extension drop lock provides 5 positions (-10°, 0°, 10°, 20°, 30°)
- Universal design fits both right and left legs
- Quick clip buckles for easy application and removal
- Comfortable padding throughout the brace

TYPICAL INDICATIONS

- ACL, PCL, MCL, LCL repairs/injuries
- Tibial plateau fractures
- Osteochondral repairs, Meniscal repairs, Patella tendon repairs
- Condylar fractures
- Sprains and strains of the knee
- HTOs (High Tibial Osteotomy)

Comfortland Universal Tri-Panel Knee Immobilizer

The Universal Tri-Panel Knee Immobilizer is a great choice for patients who require full knee immobilization. Perfect for ACL, MCL, PCL as well as soft tissue injuries, the tri-panel is fully adjustable to accommodate a wide variety of thigh and calf sizes, and is available in lengths ranging from 16"-24".

FULLY ADJUSTABLE TRI-PANEL FOR FULL KNEE IMMOBILIZATION



KEY FEATURES

- Unique design fits most thigh and calf sizes
- Sturdy construction provides full knee immobilization
- Fully adjustable tri-panel design
- Lengths range from 16" - 24"
- HCPCS L1830

TYPICAL INDICATIONS

- ACL, MCL, PCL, and LCL injuries
- Soft tissue injuries
- Post-operative knee immobilization

KNEE BRACES

Comfortland Post-Op Transition Knee Brace

Designed for post-operative knee immobilization following meniscal and chondral defect procedures, the Post-Op Transition knee brace features a telescoping technology and an independent strapping system to help provide a custom fit. Quick-release straps make it easy to don/doff, while the easy to use ROM hinge limits undesired flexion/extension. The telescoping frame along with the independent strapping system make it a universal brace, fitting patients of nearly every size.

POST-OP KNEE IMMOBILIZATION WITH ADJUSTABLE ROM HINGE



KEY FEATURES

- Universal fit
- Telescoping frame
- Quick release buckles for easy application
- Adjustable ROM hinge
- Independent strapping system
- HCPCS L1833

TYPICAL INDICATIONS

- Post-operative knee immobilization
- Tears of cartilage or meniscus of knee
- Injuries of the collateral ligaments

RCAI Knee Cage Orthosis

The Knee Cage by RCAI was designed to control Genu Recurvatum (hyperextension) of the knee. The RCAI Knee Cage permits knee flexion, aids in ambulation, and increases motion power of the quadriceps. The lightweight, thermoplastic construction of the brace conforms to the anatomy of the leg. The adjustable 3 point pressure strap behind the knee is padded for patient comfort.

LIGHTWEIGHT, THERMOPLASTIC BRACE TO CONTROL HYPEREXTENSION



KEY FEATURES

- Swedish type support
- Controls hyperextension
- Lightweight
- Padded straps for comfort
- HCPCS L1850

TYPICAL INDICATIONS

- Hyperextension of the knee due to plantar flexion contracture, quadriceps weakness, and Hip Orthoses extension weakness.

LOWER EXTREMITY



Lower extremity products for the foot and ankle.

Ankle Support

Ankle Brace

CAM walker pneumatic

CAM walker non-pneumatic

LOWER EXTREMITY

Ossur Formfit® Ankle Brace with Figure 8

The Formfit® Ankle brace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. The low-profile design fits easily into shoes and the quick lace up feature makes it easy to don and doff. Breathable material ensures it can be comfortably worn all day.

BREATHABLE MATERIALS WITH QUICK LACE UP FEATURE



KEY FEATURES

- Breathable material
- Inversion/Eversion control
- Quick lace up
- Figure 8 strap version
- HCPCS L1902

TYPICAL INDICATIONS

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

LOWER EXTREMITY

Comfortland Accord Ankle Brace III

A low-profile, lightweight ankle brace that prevents abnormal ankle inversion, eversion, and rotation while allowing natural, unrestricted dorsi and plantar-flexion. Combination of soft goods with a rigid foot plate and adjustable calf cuff provide unsurpassed levels of control and support. The Accord III also features a detachable posterior panel that stops just below the knee, giving patients the option for added support (L1971).

LOW-PROFILE, LIGHTWEIGHT ANKLE BRACE



*L1971 with panel



*L1906 without panel

KEY FEATURES

- Low profile
- Functional hinge
- Front closure design
- Rigid uprights
- Breathable, lightweight materials
- HCPCS L1971

TYPICAL INDICATIONS

- Serious sprain, injury, surgery or to stabilize, the foot or ankle fractures.
- Varus deformity, flexion deformity, osteitis deformans of ankle and foot, ankle contracture.

LOWER EXTREMITY

Ossur Formfit® Walker Air

The Formfit® Walker Air provides immobilization for patients suffering from soft tissue injury, grades 2 & 3 strains and sprains, and stable fractures. The patented pneumatic system is designed to increase stability while decreasing pain and swelling. The intuitive pump and clearly marked decompression valve make air/pressure adjustments easy.

INTEGRATED PNEUMATIC SYSTEM AND SHOCK ABSORBING SOLE



High-Top



Low-Top

KEY FEATURES

- High and Low Top versions available
- Low, wide rocker bottom promotes a stable, natural gait
- Breathable, seamless liner
- Integrated pneumatic system
- Shock absorbing sole for reduced impact at heel strike
- HCPCS Code L4361

TYPICAL INDICATIONS

- Foot and ankle conditions that may benefit from motion restriction.
- Motion restriction such as: Grade 2 and 3 ankle sprains, soft tissue injuries, stable fractures, trauma, during postoperative or rehab phase.

LOWER EXTREMITY

Formfit® Walker

The Formfit® Walker provides immobilization for patients suffering from soft tissue injury, grades 2 & 3 strains and sprains, and stable fractures.

LIGHTWEIGHT, BREATHABLE LINER WITH SHOCK ABSORBING SOLE



High-Top



Low-Top

KEY FEATURES

- High and Low Top versions available
- Shock absorbing sole for reduced impact at heel strike
- Low, wide rocker bottom promotes a stable, natural gait
- Breathable, seamless liner
- HCPCS Code L4387

TYPICAL INDICATIONS

- Foot and ankle conditions that may benefit from motion restriction.
- Motion restriction such as: Grade 2 and 3 ankle sprains, soft tissue injuries, stable fractures, trauma, during postoperative or rehab phase.

LOWER EXTREMITY

Comfortland Premium Air Walker

The premium high top and low top air walkers provide the same durability as the Comfortland CAM walkers with the addition of a pneumatic air bladder system that inflates around the ankle, providing customized compression. Designed to treat soft tissue injuries, Grade II and III sprains as well as stable fractures, the premium air walker combines functionality and superior fit into a comfortable and affordable product.

DURABLE CAM WALKER WITH PNEUMATIC AIR BLADDER SYSTEM



High-Top



Low-Top

KEY FEATURES

- Inflatable air bladder
- Shock absorbing insole
- Fits left or right foot
- HCPCS L4361

TYPICAL INDICATIONS

- Soft tissue injuries
- Grade II and III sprains
- Stable fractures

LOWER EXTREMITY

Comfortland Premium CAM Walker

The premium high top and low top CAM walker provides optimal support while incorporating many unique features that are specifically designed for both fit and comfort. The rocking walker base is designed to be low-profile and wide to promote a natural, stable gait. With its easy-to-use Velcro straps and contoured design, the cam walker minimizes impact while walking, helping to prevent further injuries.

LOW-PROFILE ROCKING WALKER BASE



High-Top



Low-Top

KEY FEATURES

- Rocker bottom promotes normal gait
- Shock absorbing insole
- Fits left or right foot
- HCPCS L4387

TYPICAL INDICATIONS

- Soft tissue injuries
- Grade II and III sprains
- Stable fractures

UPPER EXTREMITY



Upper extremity products for the wrist, elbow, and shoulder.

Wrist Brace
Elbow Brace
Shoulder Brace

Össur Formfit® Universal Wrist

The Formfit® Universal Wrist is designed to restrict motion of the wrist. The lightweight 8 inch universal design provides protection and immobilization for patients. Also available in an XL version.

UNIVERSAL FIT WITH MALLEABLE PALMAR STAY



KEY FEATURES

- Breathable materials
- Quick-pull close strap
- Positional dorsal stay
- Malleable palmar stay
- Universal 8 inch size
- Left/right specific
- HCPCS L3908

TYPICAL INDICATIONS

- Sprain and Instability
- Rheumatoid Arthritis
- Carpal Tunnel Syndrome
- Post-Op
- Fracture Stabilization
- Colles' fracture

Aspen Hinged Wrist

The Aspen Hinged Wrist is the ideal product for patients suffering from a wrist injury. The universal design allows for either right or left wrist application and the easy closure system makes it easy to don and doff. The aluminum frame can be customized to the patient's forearm, wrist and hand for a secure fit. The Aspen Hinged Wrist can be applied in either flexion or extension control to appropriately treat multiple conditions.

UNIVERSAL DESIGN FOR LEFT OR RIGHT APPLICATION



KEY FEATURES

- Universal Size
- Fits left or right
- Malleable Frame
- Adjustable wrist hinge
- HCPCS L3915/L3916

TYPICAL INDICATIONS

- Ulnar styloid fracture
- Partial Rupture of Tendon in wrist/hand
- Colles Fracture, Unspecified pain in wrist/hand
- Radius fracture/extensor tendon rupture

Comfortland Endeavor Deluxe Wrist

The Comfortland Endeavor Deluxe Wrist/Hand Splint is ideal for patients suffering from strong flexion synergy. Universally sized with an ambidextrous design, it also features a malleable frame, adjustable wrist hinge, hand bar, and finger separators to allow practitioners to achieve a perfect fit for any patient.

UNIVERSAL SIZE & AMBIDEXTEROUS DESIGN



KEY FEATURES

- Universal Size
- Fits left or right
- Malleable Frame
- Adjustable wrist hinge
- HCPCS L3915/L3916

TYPICAL INDICATIONS

- Ulnar styloid fracture
- Partial Rupture of Tendon in wrist/hand
- Colles Fracture
- Unspecified pain in wrist/hand
- Radius fracture/extensor tendon rupture

Comfortland Premium Wrist/Hand Splint

The Comfortland Premium Wrist/Hand Splint is ideal for patients suffering from carpal tunnel, osteoarthritis, rheumatoid arthritis, pain in the wrist/hand, and many other conditions. Universally sized with an ambidextrous design, it also features a malleable frame, adjustable wrist hinge, and a built-in nylon glove to aid in compression and pain relief.

UNIVERSAL DESIGN FITS LEFT & RIGHT WRIST



KEY FEATURES

- Universal Size
- Fits left or right
- Malleable Frame
- Adjustable wrist hinge
- HCPCS L3915/L3916

TYPICAL INDICATIONS

- Carpal tunnel
- Osteoarthritis
- Rheumatoid arthritis
- Pain in the wrist/hand
- Other conditions

Ossur Innovator X Post-Op Elbow Brace

With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, the Innovator X Post-Op elbow brace is designed to increase patient satisfaction and compliance.

THE MOST USER-FRIENDLY HINGE IN A POST-OP BRACE



KEY FEATURES

- Universal fit
- Left/right specific
- Comfort pad at the neck and shoulder.
- 3-dimensional padding and Lycra® lining
- OneTouch button enables easy release for slide-to-size length adjustments.
- Easy-to-use Innovator hinge makes flexion and extension adjustments simple.
- HCPCS L3760, L3761

TYPICAL INDICATIONS

- Indications requiring elbow immobilization and/or range of motion restriction may include: stable fractures of the elbow, distal humerus, proximal radius or ulna, tendon and ligament injuries or repairs, collateral ligament reconstructions, elbow hyperextension, triceps tendon repair, chronic elbow injuries, tennis elbow release.

Comfortland Premium Elbow Brace

The premium elbow brace is designed to treat any injury in which a limited range of elbow motion is needed or desired. With a lightweight design, universal sizing, and a ROM hinge with a wide variety of settings, this is a great alternative to casting and repetitive splinting. The moisture-wicking liner helps to keep patients cool and dry while the fully-adjustable design and strapping system not only make it comfortable and easy to apply, but also make the premium elbow brace one of the most cosmetically appealing elbow braces available.

PREMIUM ROM ELBOW BRACE



KEY FEATURES

- Universal size
- Left/Right specific
- Lightweight design
- Telescoping frame fits most arm lengths
- Great alternative to casting or splinting
- HCPCS L3761

TYPICAL INDICATIONS

- Cubital tunnel syndrome
- Dislocation of elbow
- Rupture of tendon in elbow

United Ortho Shoulder Abduction Sling

United Ortho's Shoulder Sling Series provides multiple options for post-operative care – immobilization, abduction, and external rotation. The innovative design creates maximum comfort and support for the patient recovering from shoulder surgery or injury. The shoulder sling is designed for a universal left or right application.

SHOULDER ABDUCTION SLING STANDARD POSITION



KEY FEATURES

- Superior fit through contoured forearm design
- Quick-release buckles for easy single handed application and removal
- Advanced moisture-wicking materials
- Premium shoulder strap padding for increased comfort
- Universal left or right
- Standard/Neutral pillow options – includes exercise ball
- HCPCS L3660

TYPICAL INDICATIONS

- Anterior glenohumeral
- Dislocation/subluxation
- Capsular shifts
- Posterior shoulder repair
- Post-operative Bankart procedures
- Soft tissue strains/repairs
- Rotator cuff repairs

Comfortland Shoulder Immobilizer w/ Waist Strap

Designed to treat rotator cuff injuries, shoulder instabilities, and anterior dislocation of the G-H joint, shoulder immobilizer features an ambidextrous design, fitting either the right or left arm. The ventilated canvas sling keeps patients cool and dry while the simple strapping design allows for easy application and adjustment.

AMBIDEXTROUS SHOULDER IMMOBILIZER



KEY FEATURES

- Ambidextrous design fits left or right arm
- Ventilated canvas construction
- Easy to apply and remove
- HCPCS L3660

TYPICAL INDICATIONS

- Rotator cuff injuries
- Shoulder instabilities
- Anterior dislocation (G-H joint)

BONE GROWTH STIMULATORS



Our bone growth stimulators treat the spine, as well as long and short bones.

**Short and Long Bone Non-Union Fracture
Upper and Lower Spine
Non-Invasive Stimulators**

BONE GROWTH STIMULATORS

DonJoy CMF Spinalogic

A portable, battery-powered, micro-controlled, noninvasive bone growth stimulator indicated as an adjunct electromagnetic treatment to primary lumbar spinal fusion surgery for one or two levels. The Spinalogic includes many features such as padded comfort, anatomical cushioning, contoured fit, and high quality materials.

NONINVASIVE ADJUNCT THERAPY BONE GROTH STIMULATOR



KEY FEATURES

- Portable, battery powered non-invasive bone growth stimulator
- Adjunct to lumbar fusion surgery
- Works with both anterior and posterior approaches
- Improves fusion rate by 49%
- Padded with contoured cushioning ensures patient comfort
- HCPCS E0748

TYPICAL INDICATIONS

- Adjunct electromagnetic treatment to primary lumbar spinal fusion surgery for one or two levels.

BONE GROWTH STIMULATORS

DonJoy CMF OL1000

For noninvasive treatment of an established nonunion fracture acquired secondary to trauma, excluding all vertebrae and flat bone. CMF OL1000 Bone Growth Stimulators are portable, battery-powered medical devices indicated for use in the noninvasive treatment of an established nonunion fracture acquired secondary to trauma, excluding all vertebrae and flat bones. Single and Dual Coil Stimulators are available in multiple sizes for fractures of short and long bones.

BONE GROWTH STIMULATOR FOR SHORT AND LONG BONES



KEY FEATURES

- For noninvasive treatment of an established nonunion fracture acquired secondary to trauma, excluding all vertebrae and flat bone.
- Lightweight and comfortable
- Easy-to-use & Noninvasive
- Requires simple, one-button operation
- Device is worn for only 30 minutes per day
- Can be used with internal/external fixation or over a cast
- HCPCS E0747

TYPICAL INDICATIONS

- Noninvasive treatment of an established nonunion acquired secondary to trauma, excluding vertebrae and flat bones. A nonunion is considered to be established when the fracture site shows no visibly progressive signs of healing.

ELITE MEDICAL SUPPLY

Elite Medical Supply dispenses quality bracing for your patients and provides superior service with online training and follow-up support to ensure patient compliance and satisfaction. We provide high-quality medical equipment, like back braces and knee braces, to patients who are being treated for cervical, spine, knee or extremity pain.

We are part of Medicare's Competitive Bidding Program for durable medical equipment, prosthetics, orthotics and supplies (DMEPOS). Under this program, Medicare chooses accredited suppliers, like us, to ensure patients have access to quality medical equipment from trusted suppliers.

Contact us today to learn more about our products and services.

www.elitemedicalsupplyofny.com

National Product Catalog v3

P 866-712-0881

To Order: order@elitemedicalsupply.com

F 866-712-0882

Questions: info@elitemedicalsupply.com

1900 Ridge Road, Suite 125

West Seneca, NY 14224